

**F**  
**Mounting holes**  
 The wall-mounted plate that holds the x-ray unit must be securely fixed to the wall. It is a suitable system for fixing the plate according to the characteristics of the wall and its ability to resist a pull force of 220 kg applied at each anchorage point.  
 If the wall is made of concrete or solid bricks, use the wall plugs supplied (m10).  
 Alternatively we recommend using the "TESCHER" chemical wall plugs which include:  
 - Broad-type injection insert (item FP 10493) - Trovabrid bar with bar and washer (item FP C 150) to

**POWER SUPPLY: 230VAC/115VAC - 50/60Hz**  
**ATTENTION: RESPECT THE CORRECT CONNECTION POLARITY.**  
**NOTE:** verify that the voltage indicated on the name plate is suitable before attempting the connection.  
**MINIMUM REQUIREMENT:**  
 Power cord 3x2.5mm<sup>2</sup> - The power cord must run at least 1 meter out of joint (K).  
 The electrical wiring must be effectively grounded in compliance with I.E.C. - us national electrical code - and C.E.I standards (or in any case all current local standards).  
 In Italy the electrical wiring must be done in compliance with electrical code C.E.I. 64-8 sez.710, that require a differential switch with the following specifications as installed before the system:  
**WITH POWER SUPPLY 230V - 50Hz**      **WITH POWER SUPPLY 115V - 60Hz**  
 NOMINAL CURRENT: 16A (RATED 0)      NOMINAL CURRENT: 16A (RATED 0)  
 DIFFERENTIAL SENSITIVITY: 0.03A      DIFFERENTIAL SENSITIVITY: 0.03A  
 A CLEARANCE OF AT LEAST 3mm

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**K**  
**POWER SUPPLY: 230VAC/115VAC - 50/60Hz**  
**ATTENTION: RESPECT THE CORRECT CONNECTION POLARITY.**  
**NOTE:** verify that the voltage indicated on the name plate is suitable before attempting the connection.  
**MINIMUM REQUIREMENT:**  
 Power cord 3x2.5mm<sup>2</sup> - The power cord must run at least 1 meter out of joint (K).  
 The electrical wiring must be effectively grounded in compliance with I.E.C. - us national electrical code - and C.E.I standards (or in any case all current local standards).  
 In Italy the electrical wiring must be done in compliance with electrical code C.E.I. 64-8 sez.710, that require a differential switch with the following specifications as installed before the system:  
**WITH POWER SUPPLY 230V - 50Hz**      **WITH POWER SUPPLY 115V - 60Hz**  
 NOMINAL CURRENT: 16A (RATED 0)      NOMINAL CURRENT: 16A (RATED 0)  
 DIFFERENTIAL SENSITIVITY: 0.03A      DIFFERENTIAL SENSITIVITY: 0.03A  
 A CLEARANCE OF AT LEAST 3mm

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